

Welcome to Laerdal Scenario Cloud



Scenario Library

Plans and pricing

Sign in

Curated scenarios.
Made by experts.
Ready for you.

Introducing **Laerdal Scenario Cloud**

[Start a free trial](#)



COVID-19 Scenarios

Train learners in how to handle a suspected case of Coronavirus (COVID-19).



<https://scenariocloud.laerdal.com/>

Getting started

Never used SimStore

SimStore Account but never used and no scenarios on account

SimStore Account and have scenario you purchased

Everyone will start by either going to our website and searching Scenario cloud or directly by going to <https://scenariocloud.laerdal.com/>

Scenario Cloud Support 1-877-523-7325



Laerdal

helping save lives

Scenario Cloud Key values

1. ***{Save time}*** *Provide access to high quality, validated content that is easy to find and run.*
2. ***{Build confidence}***
 - *Overcome the biggest technical barriers by providing a seamless workflow between Scenario Cloud and the simulation environment. (online & offline)*
 - *Encourage scenario adoption from day 1 with visual and intuitive scenario layout, supporting articles and tailored content packages.*

Laerdal Scenario Cloud Content



Find hundreds of clinical scenarios

Expand your curriculum with a selection of high-quality, expert-validated scenarios from Laerdal and partners. All available for your entire organization - with just one single license.

[Explore scenarios >](#)



Work better together

Organize scenarios in personal collections based on patient cases, curriculum goals, or other themes. In upcoming releases, you will be able to share and collaborate on collections with colleagues in your organization.



Ready when you are

Scenario allocation is now a thing of the past. Synchronize your collections across all your LLEAP and SimPad devices by simply signing in to your account.

[How to get started >](#)

Supported by our partners

American Academy
of Pediatrics
DEDICATED TO THE HEALTH OF ALL CHILDREN™



American
Heart
Association®




National League
for Nursing



Scenarios Access


- **FREE:**

- You will have access to a number of scenarios, that will typically be related to a specific simulator. We have provided these scenarios to help you understand the full value of the relevant patient simulator. This is also where we can add free scenarios when we release a new simulator or an event like the current COVID 19 where we offer free scenario. They are branded: Starter Scenarios



- **PROFESSIONAL:**

- You can get access to these highly recognized partner validated scenarios – via an institutional subscription. When your organization have a subscription, the users will have access to the entire library of scenarios.



- The scenarios can be operated simultaneously

Onboarding to Laerdal Scenario Cloud



Welcome to Laerdal Scenario Cloud

Discover what you can do with Laerdal Scenario Cloud.
It only takes 2 minutes!

Start

Overview 1. Prepare 2. **Simulate** 3. Debrief



Support during simulation

The **Simulate tab** contains all the information and files you need during a simulation session.




Next

<https://scenariocloud.laerdal.com/onboarding>

Now you are ready to start using Laerdal Scenario Cloud

Search for relevant scenario(s)



Explore Scenario



Why use this scenario?

This scenario is designed to assess a learner's ability to care for pregnant patients who present with a variety of medical conditions. They will have the opportunity to perform a point-of-care ultrasound examination using the SonoSimulator®. Learners will need to develop an appropriate differential diagnosis, initiate timely resuscitative interventions, prioritize diagnostic interventions, order appropriate basic laboratory tests, and correctly acquire, interpret, and apply point-of-care ultrasound findings towards medical decision-making.

Move to collection(s)

My Scenarios

Collections

My Default Collection
3 scenarios

My Previous Purchased
2 scenarios

Covid 19
6 scenarios

New collection

My Default Collection



Resuscitation of 25-Week Gestation Newborn with PPV, CPAP, Intubation, and Surfactant Administration

AMERICAN ACADEMY OF PEDIATRICS



Sepsis secondary to urinary tract infection

LAERDAL MEDICAL



Obstructive Uropathy from Ureteral Calculus, Intrauterine Pregnancy of Uncertain Viability

SONOSIM


Synchronize automatically My collection – LLEAP - SimPad

Sign on to Laerdal Scenario Cloud from LLEAP - SimPad




How to operate from LLEAP - SimPad

How to operate. Software version 7.0 or higher



How to update your devices

Play the videos below to learn how to update SimPad and/or LLEAP to the latest software version.



[View instructions for updating SimPad](#)



[View instructions for updating LLEAP](#)

Scenario Cloud Help

Update your equipment

To be able to use scenarios you must **update to the latest version of LLEAP and SimPad**. We have created videos to help guide you through this process.

[Updating your device \(video guide\)](#)

How to access scenarios on SimPad

Sign in to your Laerdal account on SimPad to access the collections you create on Laerdal Scenario Cloud.

[Accessing scenario collections on SimPad](#)



How to access scenarios using LLEAP

Sign in to your Laerdal account in LLEAP to access your collections and run scenarios on LLEAP.


[Accessing scenario collections on LLEAP](#)

<https://scenariocloud.laerdal.com/help>

Sign on to Laerdal Scenario Cloud from LLEAP - SimPad



Seamless flow from Scenario Cloud to operation




Sign on to Laerdal Scenario Cloud from LLEAP - SimPad



- **Log on to Laerdal Scenario Cloud**

- Transfer selected scenarios from general library into my scenarios
- Creating your personal collection(s)



- **Log on to Laerdal Scenario Cloud**

- Use your personal log – in from your SimPad and LLEAP devices
- While signed in the collections from Scenario Library will be synchronized automatically
- You can now run the scenarios

CARES Act

- **Cover Letter from CARES**
<https://documentcloud.adobe.com/link/track?uri=urn%3Aaaid%3AscDs%3AUS%3A6390bc57-418e-499d-9bd1-997e8b495777>
- **Frequently Asked Questions**
<https://documentcloud.adobe.com/link/track?uri=urn%3Aaaid%3AscDs%3AUS%3A0876bcd4-8dc7-4c84-a44c-d93f0aa63cb4>
- **Methodology**
<https://documentcloud.adobe.com/link/track?uri=urn%3Aaaid%3AscDs%3AUS%3A44866b6e-62dc-46d9-8866-6bb308943d79>
- **Agreement**
<https://documentcloud.adobe.com/link/track?uri=urn%3Aaaid%3AscDs%3AUS%3A12bce09f-ef9c-476f-a73c-3442fdffe495>
- **Allocation**
<https://documentcloud.adobe.com/link/track?uri=urn%3Aaaid%3AscDs%3AUS%3A4bdec716-6103-4d57-953c-24bead080295>
- Thanks



Laerdal

helping save lives