


## Downloading the vrClinicals for Nursing App for Meta Quest 2 Headset


You can download the app directly from the headset, a browser, or the Meta Quest mobile app

### Download from Headset

1. Start your Meta Quest 2 headset and put it on your head
2. From the main menu, select the App Library icon



3. Select the Get More Apps option





Wolters Kluwer




Laerdal


helping save lives

4. Type in "vrclinicals" in the search bar at the top and select the Enter key


**NOTE: You must type in the exact spelling, do not put a space between vr and clinicals**



5. Select the vrClinicals app icon



6. Select View Details






Wolters Kluwer



Laerdal

helping save lives


7. Select OK on the App Lab notice



8. Select Download



9. Return to the App Library






Wolters Kluwer



Laerdal

helping save lives

10. The app will be displayed and download progress will be indicated



11. Once the download is complete, the app is available in your app library






Wolters Kluwer




Laerdal  
helping save lives

## Download from Browser


1. On your computer, navigate to [www.oculus.com/experiences/quest/](https://www.oculus.com/experiences/quest/)
2. Log in to the same Meta account that your desired headset is logged in to



3. Type “vrClinicals” in the search bar and hit Enter



4. On the search results page, click on View App






Wolters Kluwer




Laerdal

helping save lives


- Click on the icon for the app



- Click on OK on the App Lab message



- Click on the Get button



- Put your headset on, navigate to the app library, and install the vrClinicals for Nursing app




Wolters Kluwer



Laerdal  
helping save lives

## Download from Mobile Device

1. Open the Meta Quest app on your mobile device
2. Ensure you are logged in to the same Meta account that your desired headset is logged in to
3. Select the Search icon from the bottom menu





Wolters Kluwer




Laerdal


helping save lives

4. Type in “vrclinicals” in the search field and click search


**NOTE: You must type in the exact spelling, do not put a space between vr and clinicals**




5. Press on the vrClinicals app icon



6. Select View App from the App Lab message



7. Select OK on the pop-up modal and then select Get



8. Put your headset on, navigate to the app library, and install the vrClinicals for Nursing app